

**NAME**

pyFAI-recalib – pyFAI-recalib

**SYNOPSIS****pyFAI-recalib** [*options*] -p *ponifile* -w *l* -S *calibrant.D* *imagefile.edf***DESCRIPTION**

Calibrate the diffraction setup geometry based on Debye–Sherrer rings images with a priori knowledge of your setup (an input PONI–file). You will need a "d–spacing" file containing the spacing of Miller plans in Angstrom (in decreasing order). If you are using a standart calibrant, look at <https://github.com/kif/pyFAI/tree/master/calibration> or search in the American Mineralogist database: <http://rruff.geo.arizona.edu/AMS/amcsd.php>

**OPTIONS****--version**

show program's version number and exit

**-h, --help**

show this help message and exit

**-o FILE, --out=FILE**

Filename where processed image is saved

**-v, --verbose**

switch to debug/verbose mode

**-S FILE, --spacing=FILE**

file containing d–spacing of the reference sample (MANDATORY)

**-w WAVELENGTH, --wavelength=WAVELENGTH**

wavelength of the X–Ray beam in Angstrom

**-e ENERGY, --energy=ENERGY**

energy of the X–Ray beam in keV (hc=12.398419292keV.A)

**-P POLARIZATION\_FACTOR, --polarization=POLARIZATION\_FACTOR**polarization factor, from **-1** (vertical) to **+1** (horizontal), default is None (no correction), synchrotrons are around 0.95**-b BACKGROUND, --background=BACKGROUND**

Automatic background subtraction if no value are provided

**-d DARK, --dark=DARK**

list of dark images to average and subtract

**-f FLAT, --flat=FLAT**

list of flat images to average and divide

**-s SPLINE, --spline=SPLINE**

spline file describing the detector distortion

**-D DETECTOR\_NAME, --detector=DETECTOR\_NAME**

Detector name (instead of pixel size+spline)

**-m MASK, --mask=MASK**

file containing the mask (for image reconstruction)

**-n NPT, --pt=NPT**

file with datapoints saved. Default: basename.npt

**--filter=FILTER**

select the filter, either mean(default), max or median

**-l DISTANCE, --distance=DISTANCE**

sample–detector distance in millimeter

**--poni1=***PONI1*  
poni1 coordinate in meter

**--poni2=***PONI2*  
poni2 coordinate in meter

**--rot1=***ROT1*  
rot1 in radians

**--rot2=***ROT2*  
rot2 in radians

**--rot3=***ROT3*  
rot3 in radians

**--fix-dist**  
fix the distance parameter

**--free-dist**  
free the distance parameter

**--fix-poni1**  
fix the poni1 parameter

**--free-poni1**  
free the poni1 parameter

**--fix-poni2**  
fix the poni2 parameter

**--free-poni2**  
free the poni2 parameter

**--fix-rot1**  
fix the rot1 parameter

**--free-rot1**  
free the rot1 parameter

**--fix-rot2**  
fix the rot2 parameter

**--free-rot2**  
free the rot2 parameter

**--fix-rot3**  
fix the rot3 parameter

**--free-rot3**  
free the rot3 parameter

**--fix-wavelength**  
fix the wavelength parameter

**--free-wavelength**  
free the wavelength parameter

**--saturation=***SATURATION*  
consider all pixel>max\*(1-saturation) as saturated and reconstruct them

**--weighted**  
weight fit by intensity, by default not.

**--npt=***NPT\_1D*  
Number of point in 1D integrated pattern, Default: 1024

- npt-azim**=*NPT\_2D\_AZIM*  
Number of azimuthal sectors in 2D integrated images. Default: 360
- npt-rad**=*NPT\_2D\_RAD*  
Number of radial bins in 2D integrated images. Default: 400
- unit**=*UNIT*  
Valid units for radial range: 2th\_deg, 2th\_rad, q\_nm<sup>-1</sup>, q\_A<sup>-1</sup>, r\_mm. Default: 2th\_deg
- no-gui**  
force the program to run without a Graphical interface
- no-interactive**  
force the program to run and exit without prompting for refinements
- r** *MAX\_RINGS*, **--ring**=*MAX\_RINGS*  
maximum number of rings to extract. Default: all accessible
- p** *FILE*, **--poni**=*FILE*  
file containing the diffraction parameter (poni-file). MANDATORY
- k**, **--keep**  
Keep existing control point and append new

The main difference with pyFAI-calib is the way control-point hence DebyeSherrer rings are extracted. While pyFAI-calib relies on the contiguity of a region of peaks called massif; pyFAI-recalib knows approximately the geometry and is able to select the region where the ring should be. From this region it selects automatically the various peaks; making pyFAI-recalib able to run without graphical interface and without human intervention (**--no-gui** **--nointeractive** options).